


Specific Instructions for executing the RNASeq pipeline

Introduction


The RNASeq pipeline has been implemented to quantify gene & isoform expression and carry out tests for differential expression at both the gene- and isoform- level. The pipeline also generates results for differential expression across splice-junctions as well as differential expression at the gene-level, counting only junction-spanning reads.



The RNASeq pipeline is run in two phases. First all reads must pass through the initialQC phase of adapter-clipping, read-level QC, read-alignment and alignment-level. Second genes and transcripts counts and differential expression analysis are generated.

Running the RNASeq pipeline

Please refer to the CCBR Pipeliner overview documentation for specific instructions to launch the CCBR pipeliner. Here we assume you succeed to launch the GUI as shown below



The screenshot shows a window titled "CCBR Pipeliner" with a standard macOS-style title bar (red, yellow, green buttons). The window has a menu bar with "File", "View", and "Help". The main content area is divided into three sections:

- Project Information:** This section contains three input fields with labels and examples:
 - Project Id:** The input field contains "project". To its right, the text "(Examples: CCBR-nnn, Labname or short project name)" is displayed.
 - Email address:** The input field is empty. To its right, the text "(Mandatory field: must use @nih.gov email address)" is displayed.
 - Flow Cell ID:** The input field contains "stats". To its right, the text "(Examples: FlowCellID, Labname, date or short project name)" is displayed.
- Global Settings:** This section contains two dropdown menus:
 - Genome:** The dropdown menu is set to "hg19".
 - Pipeline Family:** The dropdown menu is set to "Select a pipeline".
- Project Description:** This section features a large text area for input. The text "Enter CCBR Project Description and Notes here." is at the top. On the right side of the text area, there are three small, vertically stacked buttons (up, down, and a middle button).

Phase1: Initial QC

After entering information on the 'Project Information and 'Global Options' tab, select the 'rnaseq' pipeline from the pipeline family list box and set the RNASeq Pipeline Options as follows:

- 1) First, select or type the full path for the data and working directory. The working directory name **should be new** and the **system will create it automatically**
- 2) Select the sub-pipeline “**Quality Control Analysis**” and ignore the sample information that are not used for the Initial QC phase.
- 3) Click on button “initialize Directory” to initialize the working directory
- 4) Click on button “Dry run” to check the initialqc phase tasks and settings.
- 5) **If there is no error message in the “Dry run”**, then click on the button “Run” to run this pipeline QC phase on Biowulf cluster.

NOTE: Once the pipeline job is submitted to Biowulf cluster, you will be notified by email and you can close the GUI. On the command prompt, you can check the status of your jobs by using this command:

➤ `sjobs (or squeue -u <username>)`

CCBR Pipeliner

File View Help

Project Information

Project Id: maseq_01 (Examples: CCB-NNN, Labname or short project name)

Email address: user@mail.nih.gov (Mandatory field: must use @nih.gov email address)

Flow Cell ID: FlowcellIXY (Examples: FlowCellID, Labname, date or short project name)

Global Settings

Genome: hg19 Pipeline Family: maseq

Project Description RNAseq 1

Data Directory: [] Open Directory

FastQ files Found: 0

Working Directory: [] Open Directory

3 Initialize Directory 4 Dry Run 5 Run

Options 2

Pipeline: Quality Control Analysis

Sample Information

Set Groups Set Contrasts

Results from Initial QC

Once the pipeline completes successfully, the working directories will contain the following folders/files:

- **QC folder:** contains Fastqc results on the adapter-clipped fastq files
- **trim folder:** contains temporary adapter-clipped fastq files (both paired and unpaired)
- **STAR 1st pass output files:**
 - *.Aligned.out.bam
 - *.SJ.out.tab
 - *.Log.progress.out
 - *.Log.final.out
- **STAR 2nd pass output files:**
 - *.p2.Aligned.toTranscriptome.out.bam
 - *.p2.ReadsPerGene.out.tab
 - *.p2.Aligned.sortedByCoord.out.bam
 - *.p2.SJ.out.tab
 - *.p2.Log.progress.out
 - *.p2.Log.final.out
 - *.star_rg_added.sorted.bam
 - *.star_rg_added.sorted.dmark.bam
 - *.star_rg_added.sorted.dmark.bai
- **Alignment QC files:**
 - a) *.RnaSeqMetrics.txt: Output from Picard CollectRNASeq Metrics
 - b) *.flagstat.concord.txt: Output from samtools flagstat
 - c) RSEQC output:
 - *.inner_distance_plot.pdf
 - *.inner_distance_plot.r
 - *.inner_distance_freq.txt
 - *.inner_distance.txt
 - *.GC.xls
 - *.GC_plot.pdf
 - *.strand.info
 - *.Rdist.info
 - d) ProjectID_FlowCellID_Summary.txt : comprehensive Alignment QC metrics from RSEQC
 - e) ProjectID_FlowCellID.xlsx: QC-metrics and barcharts in excel format

- f) FQscreen folder: contains results from Fastq-Screen
- **Reports folder:**
 - a) `aggregate_fastqc_report.html`: quick view of PASS/WARN/FAIL flags using color codes, across all samples and all metrics
 - b) **`multiqc_report.html`**: comprehensive and interactive HTML page aggregating QC metrics across several tools

Phase2: Differential expression workflow (after Initial QC)

This workflow is to be executed in the same working directory as that specified in the InitialQC phase, so that the pipeline can identify the output files from the InitialQC run, that are being used as input files to the differential expression workflow.

After entering information on the 'Project Information' and 'Global Options' tab, select the 'Differential Expression Analysis' pipeline from the pipeline family list box and set the RNASeq Pipeline Options as follows:

- 1) First, select or type the full path for the data and working directory. The working directory should be the same one you specified in the initial QC phase
- 2) Select the sub-pipeline "**Differential Expression Analysis**"
- 3) choose option "yes, Report Differentially Expressed genes" if you have a set of contrasts or choose the default "no, Do not Report Differentially Expressed genes" if only interested in genes or transcripts counts
- 4) set your criteria for gene filtering: minimum reads counts shared by at least a number of samples. The minimum read count is used to compute a minimum count per million (cpm) threshold based on the maximum library size. This filtering removes genes that are not expressed in any sample (all counts are zeros), as well as filters out genes with very low cpm below the computed cpm threshold.
- 5) click on the button "Set Groups" to define the list of samples to include in this phase with the the group (condition) and label (for plots) information. You can type all information needed. After filling in the information, make sure to click on 'Save' to save this information (a file `groups.tab` will be created). You can also load (read) the data from a predefined file "`groups.tab`". **No field name is required.**

Here one example:

<sample>	<condition>	<label>
SS.SRR950078	G1	s78
SS.SRR950079	G1	s79

SS.SRR950080	G2	s80
SS.SRR950081	G2	s81

6) click on the button “Set Contrasts” if only you have to report Report differentially expressed genes. You need to enter the groups to compare (contrasts). You can type a contrast per line where the first group (group1) will be compared to the second group (group2 as control). Make sure to save this information (a file contrasts.tab will be created). You can also load (read) the data from a predefined file “contrasts.tab”. No field name is required.

Here one example:

G2 G1

7)Click on button “Dry run” to check this phase tasks and settings

8) **If there is no error message in the “Dry run”**, then click on the button “Run” to run this pipeline differential expression analysis phase on Biowulf cluster.

CCBR Pipeliner

File View Help

Project Information

Project Id: (Examples: CCBR-nnn, Labname or short project name)

Email address: (Mandatory field: must use @nih.gov email address)

Flow Cell ID: (Examples: FlowCellID, Labname, date or short project name)

Global Settings

Genome: Pipeline Family:

Project Description **RNAseq** 1

Data Directory:

FastQ files Found: 0

Working Directory:

7 8

Options

Pipeline 2

3

4 Low Abundance Gene Thresholds

Include genes with \geq read counts in \geq samples

5-6 Sample Information

Here is a snapshot of how the “Groups” and “Contrasts” text-box would look after filling in the information:

The image displays two screenshots of the CCBR Pipeliner software interface, showing the 'Groups Information' and 'Contrasts Information' windows.

Groups Information Window:

SS	Group	Contrast
SRR950078	G1	s78
SRR950079	G1	s79
SRR950080	G2	s80
SRR950081	G2	s81S

Buttons: Load, Save

Contrasts Information Window:

Contrast	Group
G2	G1

Buttons: Load, Save

Results from Differential expression analysis

Once the pipeline completes successfully, the working directory should contain the following files/folders:

- **Folder DEG_genes:** contains results of differential expression at the gene-level based on FeatureCounts method raw counts (Subread) from three DE algorithms: Limma, DESeq2 and EdgeR
- **Folder DEG_rsemgenes:** contains results of differential expression at the gene-level based on RSEM raw counts from three DE algorithms: Limma, DESeq2 and EdgeR
- **Folder DEG_geneJunctions:** contains results of differential expression at the gene level, counting only exon-exon junction-spanning reads (useful if the library is total RNA, where we see a large proportion of intronic reads, originating from unspliced RNA). All 3 DE methods are applied to find differentially expressed genes.
- **Folder DEG_junctions:** contains results of differential expression at the junction-level. Here we are simply quantifying reads across each known splice-junction in the chosen genome. So the output contains read counts for each exon-exon junction in the genome (not collapsing by gene). Then all 3 methods are applied to find differentially expressed junctions.
- **RSEM results:**
 - *.rsem.genes.results: gene-level counts from RSEM
 - *.rsem.isoforms.results: isoform-level counts from RSEM
- **EBSeq results:**
 - *.ebseq: differential expression results from EBSeq at the isoform-level (from RSEM counts)
- **Folder salmonrun:** Contains one folder for each sample, containing salmon quantification files as well as differential expression results from Sleuth at the isoform-level